

# SMART LIBRARIES



## INSTRUCTIONS

1. DOWNLOAD SCANBODIES LIBRARIES INTO EXOCAD SOFTWARE


---

2. INSTALL SCANBODIES LIBRARIES

---


3. INSTRUCTIONS FOR Ti BASE AND Co-Cr BASE WORK

## 1.DOWNLOAD SCANBODIES LIBRARIES INTO EXOCAD SOFTWARE




### STEP 1

1. From a computer, connect to: [www.smartimplantsolutions.com](http://www.smartimplantsolutions.com)
2. Click the **LIBRARIES** button in the main menu.



### STEP 2

1. Click on the library button to download.  
This link takes you directly to the **DROPBOX** download page.



### STEP 3

#### DROPBOX:


1. Press the **DOWNLOAD** button  
The file will be downloaded to the Downloads folder on your PC.

## 2.INSTALL SCANBODIES LIBRARIES

### STEP 1

#### UNZIP THE FILE .ZIP

1. Go to the Downloads folder on your PC.
2. Unzip the downloaded file.




To unzip the file you can use a standard zip file decompression program downloaded from the internet such as WinZip or WinRar.

### STEP 2

When you unzip the file, you will see a folder called "EXOCAD" that contains:



EXOCAD



#### **IMPORTANT\_ Read before using the libraries.txt**

When you download the libraries for the first time, it is recommended do a test before using the libraries in definitive works.

.....  
Before using the special XS and XXS libraries, please, contact our technical support.

### SCANBODY Ti LIBRARIES



#### SCANBODY Ti

Libraries for Ti-Base, CrCo base (overcastable) and digital analog, from intraoral or desktop scanning.

#### SCANBODY D (Implant Level)

Libraries for titanium scanbodies. They are used for implant direct work over digital analog, from intraoral or desktop scanning.

**IMPORTANT:** The orientation of the flat face must always be towards vestibular side.



#### SCANBODY INY

Libraries to work with free scanbody that comes with every Ti-Base.

It allows you to work with Ti-Base or CrCo base.

You have to scan always with standard height Ti-Base (the lowest one).

If you want to choose different gingival or cementation heights, you have to select them inside exocad program.

## 2.INSTALL SCANBODIES LIBRARIES


### STEP 3

#### COPY FOLDERS

Copy the folders or folders from the Libraries that interest you.

Example:

1. Open the folder **SIS-IMPLANT LIBRARY\_IMPLANT LEVEL\_Intraoral\_TI....**
2. Select and copy all folders.
3. Paste them into the following patha **C:/exocad/DentalCADApp/library/implant/**




**SHORTCUT:** Using file manager of your computer, you can find required path very quickly.

**C:/exocad/DentalCADApp/library/implant/**

### STEP 4 CHECK

1. Make sure you have the libraries correctly copied in the Exocad library menu.
2. Go into **implant** folder and be sure that files are copied

### STEP 5 EXECUTE EXOCAD



1. **Execute** Exocad.
2. **Choose** the SIS Ti-Base.

### 3. INSTRUCTIONS FOR Ti BASE AND Co-Cr BASE WORK

#### STEP 1

Location of implants and connections in libraries.

#### PRODUCT LABEL



#### PRODUCT REFERENCE

**SIS-1100: Straumann® Bone Level®**

**NOTA:** The location of libraries corresponds to the same reference on the blister in hundreds.

- SIS-0100: Nobel Biocare® Branemark®
- SIS-0120: BTI® Hexágono externo
- SIS-0200: Biomet 3i® Osseotite®
- SIS-0300: Zimmer® Screw-vent®
- SIS-0400: Biomet 3i® Certain®
- SIS-0500: Astra® Tech Osseospeed®
- SIS-0600: Straumann® Tissue Level®
- SIS-0700: BTI® Tetralobular
- SIS-0800: Nobel Biocare® Replace Select®
- SIS-0900: Nobel Biocare® Nobelactive®
- SIS-1000: Klockner® Sk2®
- SIS-1100: Straumann® Bone Level®
- SIS-1300: Multi-Unit
- SIS-1400: Sweden & Martina® Premium Kohno®

- SIS-1500: Phibo® Tsh®
- SIS-1600: Biohorizons®
- SIS-1800: Mis® C1/V3
- SIS-1900: Klockner® EC
- SIS-2200: Neobiotech® Is System
- SIS-2300: Astra® Tech System Evolution
- SIS-2400: Neodent GM®
- SIS-2500: Megagen Anyridge®
- SIS-2600: Denstply® Xive®
- SIS-2700: Denstply® Ankylos®
- SIS-2800: Osstem® TS/ Hiossen® ET
- SIS-2900: Dio® UFII®
- SIS-3000: Dentium Superline & Implantium®
- SIS-3200: Klockner® Vega®

### 3. INSTRUCTIONS FOR Ti BASE AND Co-Cr BASE WORK

---

#### STEP 2

Once library is chosen, **select the desired platform.**

---

#### STEP 3

Once library is chosen, **select the desired platform.**

Later, **choose** if you want to work **over Ti-base or CrCo base.**

If you **select Ti-base option** and you want to use a **different height, cement o gingival, choose** desired option:



**HG2**

Gingival height 2mm.



**HG3**

Gingival height 3 mm.



**HC7**

Gingival height 7mm.

**IMPORTANT:** These Ti-Bases are sold without Scanbody, as they don't have a Scanbody function.

---

#### STEP 4

When you choose the abutment, you can select option:

E (Engaging) or NE (Non-Engaging )

---

#### STEP 5

In any option you can select **different frictions** (space for cement): **S, N, B, XS** or **XXS**.

**XXS: Extra Extra Small**, (less space for cement than XS, more friction).

**XS: Extra Small**, (less space for cement than S, more friction).

**S: Small** (less space for cement, more friction).

**N: Normal** (standard space for cement).

**B: Big** (more space for cement, less friction).

**Select the option you prefer.**

**IMPORTANT:** It's recommended to do a friction test before carrying out the final work.